

Product Datasheet

Name: Mab-C-P

Description: Hybridoma clones have been derived from hybridization of myeloma cells with spleen cells of BALB/c mouse immunized with purified C-Peptide.

Catalog No.	Isotype	pI	Unit	Buffer	Usage
M1021	IgG1	7.3	mg	0.01M PB+0.1M NaCl (pH7.4)	Detection

Specificity: Mabs react with C-Peptide

Format: Liquid

Concentration: ≥ 2 mg/ml

Purity: $\geq 90\%$ (SDS-PAGE)

Preservative: 0.1% proclin300

Application: Recommended for sandwich immunoassays in ELISA, CLIA, GICA. Each laboratory should determine an optimum working titer for use in its particular application.

Storage: Store at 2-8°C up to 21 days, -20°C for long term. Avoid multiple freeze-thaw cycles.

Note: Test systems and platform change may cause different results, please contact us if you have any questions.

For research or further manufacturing use only.

Nanjing Genstars Biotech Co.,Ltd.

Address: Layer 4, building A6&A7 Maple Science Park, Qixia District, Nanjing, China

Tel: +86-025-86793889

Email: marketing@genstars.com

Website: www.genstars.com



Website



Wechat